

# ALLERGENE

Schnelltest – Testdauer 10 Minuten

**AlerTox**



**AlerTox** Teststreifen für Lebensmittel und Oberflächen

Kit	Reference	LOD	Format
AlerTox Sticks Casein	KT-5781	2.5 ppm	10 strips
AlerTox Sticks Casein	KT-5772	2.5 ppm	25 strips
AlerTox Sticks $\beta$ -Lactoglobulin	KT-5782	2.5 ppm	10 strips
AlerTox Sticks $\beta$ -Lactoglobulin	KT-5773	2.5 ppm	25 strips
AlerTox Sticks Crustacea	KT-6126	7 ppm (dry) / 33 ppm (wet)	10 strips
AlerTox Sticks Peanut	KT-6128	1 ppm	10 strips
AlerTox Sticks Hazelnut	KT-6129	1 ppm	10 strips
AlerTox Sticks Almond	KT-6127	10 ppm	10 strips
AlerTox Sticks Soy	KT-6125	1 ppm	10 strips
AlerTox Sticks Egg	KT-5899	1 ppm	10 strips
AlerTox Sticks Egg	KT-5898	1 ppm	25 strips
AlerTox Sticks Fish	KT-6131	0.35 ppm (dry) / 1 ppm (wet)	10 strips
AlerTox Sticks Mustard Seeds	KT-6370	2 ppm	10 strips