

Buffer Tablets

Loba Chemie pH Buffer Tablets offers easier way to prepare buffer solutions of required concentrations. One has to simply dissolve tablet to 100 ml deionized or distilled water and a solution of the specified pH will be ready for use.

Beside buffer tablets, Loba Chemie product array also includes ready to use buffer solutions and concentrate buffer solutions with different pH.



1 to 11 of 11 item(s) displayed

50 per page

Code	Product Name	UN No.	IMCO	CAS No.	MSDS	COAs
02374	BUFFER TABLETS pH 4.0 (1 TABLET DISSOLVED TO 100 ML SOLUTION)	-	-	NA	MSDS	COAs
2374D	BUFFER TABLETS pH 6.0 (1 TABLET DISSOLVED TO 100 ML SOLUTION)	-	-	NA	MSDS	COAs
2374G	BUFFER TABLETS pH 6.40 (1 TABLET DISSOLVED TO 100 ML SOLUTION) NEW	-	-			COAs
2374H	BUFFER TABLETS pH 6.40 (1 TABLET DISSOLVED TO 1000 ML SOLUTION) NEW	-	-			COAs
2374E	BUFFER TABLETS pH 6.8 (1 TABLET DISSOLVED TO 100 ML SOLUTION)	-	-	NA	MSDS	COAs
2374F	BUFFER TABLETS PH 6.80 (1 TABLET DISSOLVED TO 1000 ML SOLUTION)	-	-			COAs
02375	BUFFER TABLETS pH 7.0 (1 TABLET DISSOLVED TO 100 ML SOLUTION)	-	-	NA	MSDS	COAs
2375A	BUFFER TABLETS pH 7.0 (10 TABLETS DISSOLVED TO 100 ML SOLUTION)	-	-	NA	MSDS	COAs
2375B	BUFFER TABLETS PH 7.00 (1 TABLET DISSOLVED TO 1000 ML SOLUTION)	-	-			COAs
2375E	BUFFER TABLETS PH 7.20 (1 TABLET DISSOLVED TO 1000 ML SOLUTION)	-	-			COAs
02376	BUFFER TABLETS pH 9.2 (1 TABLET DISSOLVED TO 100 ML SOLUTION)	-	-	NA	MSDS	COAs